


Polyester Tipped Applicator

Aluminum/plastic handle



Description

This specialty applicator is tipped with medical grade polyester that we card ourselves. They're non-absorbent and contain no fatty acids or latent sugars, making them ideal for transport, lab work, and diagnostic kits. The miniature polyester tip is adhered to flexible aluminum which is inserted into a rigid plastic handle.

Dimensions

6" Overall Length (152.4mm)

Tip

■ 0.078" Dia. (1.981mm)

■ 0.688" Length (17.475mm)

Shaft

■ 0.098" / 0.035" Dia. (2.489mm / 0.889mm)

■ 4" / 2" Length (101.6mm / 50.8mm)

Standards

- CE marked for medical use
- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of sterility and/or compliance available with each shipment upon request

Qualities

- Medical grade quality
- Synthetic polyester contains no fatty acids or latent sugars
- Thin flexible aluminum on rigid plastic handle

Applications

- Specimen collection
- Testing and screening swab

Ref #	Packaging
800-PD ALUM <i>Non-sterile</i>	10/100 (10 bags per case, 100 per bag)
25-800 1PD ALUM 50 <i>Sterile</i>	10/50/1 (10 boxes of 50 individually wrapped packages; 1 per package)



Quality Medical Products Since 1919

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

1-800-321-2313
puritanmedproducts.com