


Polyester Tipped Applicator

Aluminum shaft



- Polyester tip

- Flexible aluminum handle

Actual Size

Description

This specialty applicator is tipped with medical grade polyester that we card ourselves. They're non-absorbent and contain no fatty acids or latent sugars, making them ideal for transport, lab work, and diagnostic kits.

Dimensions

5.437" Overall Length (138.1mm)

Tip

- 0.078" Dia. (1.981mm)
- 0.687" Length (17.45mm)

Shaft

- 0.035" Dia. (0.89mm)
- 5.375" Length (136.525mm)

Standards

- CE marked for medical use
- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of sterility and/or compliance available with each shipment upon request

Qualities

- Medical grade quality
- Synthetic polyester contains no fatty acids or latent sugars
- Thin flexible aluminum handle follows contours of passages with virtually no injury to mucosal tissues

Applications

- Forensic evidence collection
- Diagnostic specimen collection for testing and screening
- Flu and strep testing
- Pediatric
- Eye microbiology
- Veterinary care

Ref #	Packaging
25-800 D 50 (sterile)	10/50/1 (10 boxes of 50 individually wrapped packages; 1 per package)



Quality Medical Products Since 1919

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

1-800-321-2313
puritanmedproducts.com