

Rayon Tipped Applicator

Aluminum handle



Actual Size

Description

This tipped applicator is made of the finest USP grade of rayon fiber that is carded in our factory to maintain consistency. Ideal for diagnostic specimen collection, wound treatment, and application of medicines. Flexible aluminum handle for minimal patient discomfort.

Dimensions

5.437" Overall Length (138.1mm)

Tip

- 0.078" Dia. (1.98mm)
- 0.687" Length (17.45mm)

Shaft

- 0.035" Dia. (0.89mm)
- 5.375" Length (136.53mm)

Standards

- CE marked for medical use
- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of compliance available with each shipment upon request

Qualities

- Medical grade quality
- Flexible aluminum handle

Applications

- Specimen collection
- General patient care
- Wound treatment
- Applying medicines

Ref #	Packaging
25-800 R 50 <i>Sterile</i>	10/50/1 (10 boxes per case, 50 individually wrapped packages, 1 per package)



Quality Medical Products Since 1919

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

1-800-321-2313
puritanmedproducts.com