


Calgiswab® Calcium Alginate Tipped Applicator

Aluminum shaft



- Calcium alginate tip

- Flexible aluminum handle

Actual Size

Description

This specialty applicator is tipped with calcium alginate fiber, an organic material derived from seaweed and carded in Puritan's facility to assure conformance with medical device specifications. Best known as an Urethro-genital swab.

Dimensions

5.437" Overall Length (138.1mm)

Tip

- 0.078" Dia. (1.981mm)
- 0.437" Length (11.1mm)

Shaft

- 0.025" Dia. (0.635mm)
- 5.375" Length (136.525mm)

Standards

- CE marked for medical use
- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of sterility and/or compliance available with each shipment upon request

Qualities

- Medical grade quality
- Calcium alginate tip contains no fatty acids, which may be inhibitory to many strains of bacteria
- Thin flexible aluminum handle for minimal patient discomfort

Applications

- Flexible urethro-genital swab
- Specimen collection for testing and screening
- STD testing
- Pediatrics
- Ophthalmology
- Veterinary care

Ref #	Packaging
25-801 A 50 (sterile)	10/50/1 (10 boxes of 50 individually wrapped packages; 1 per package)



Quality Medical Products Since 1919

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

1-800-321-2313
puritanmedproducts.com