


Polyester Tipped Applicator

Blue polystyrene handle



Actual Size

Description

This synthetic tipped applicator is made of medical-grade polyester fiber that is carded in our factory to maintain consistency. It contains no fatty acids or latent sugars, making it ideal for specimen transport, lab work, and diagnostic kits. It is bonded securely to a plastic handle by an aqueous based adhesive. Colored handle for color coded swab uses.

Dimensions

3" Overall Length (76.2 mm)

Tip

■ 0.203" Dia. (5.156 mm)

■ 0.687" Length (17.45 mm)

Handle

■ 0.100" Dia. (2.54 mm)

■ 2.937" Length (74.6 mm)

Standards

- CE marked for medical use
- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of sterility and/or compliance available with each shipment upon request

Qualities

- Medical grade quality
- Absorbent tip
- Soft smooth tip
- Colored handle for color coded uses

Applications

Medical / Diagnostic

- Paternity testing
- DNA collection
- Diagnostic testing

Ref #	Packaging
25-803 2PD BLUE <i>Sterile</i>	10/100/2 (10 boxes of 100 wrapped packages; 2 per package)



Quality Medical Products Since 1919

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

1-800-321-2313
puritanmedproducts.com