


Fab-Swab Cotton Tipped Applicator

Polystyrene handle



Not actual size

Description

This DNA controlled absorbent cotton tipped applicator is produced with High Quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a polystyrene handle.

The convenient, reclosable, dry transport device is designed for reliable sample and evidence collection. The patented tube design entails a breathable filter label that covers holes to prohibit mold growth during transport. ID labels are included.

Dimensions

5" Overall Length (127 mm)

Tip

- 0.203" Dia. (5.156 mm)
- 0.687" Length (17.45 mm)

Handle

- 0.100" Dia. (2.54 mm)
- 4.875" Length (123.825 mm)

Standards

- CE marked for medical use
- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of compliance available with each shipment upon request

Qualities

- Medical grade quality
- Produced in a DNA controlled environment
- Reclosable transport device
- Can be used dry, with IPA, or sterile water
- Includes evidence information label

Applications

- DNA testing
- Forensic evidence collection

| Ref # | Packaging |
|---|---|
| 25-805 1PC TT FABUSA <i>Sterile</i> | 10/50/1 (10 boxes per case, 50 individually wrapped tubes per box, 1 swab per box, 1 swab per tube) |



Quality Medical Products Since 1919

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

1-800-321-2313
puritanmedproducts.com