


Polyester Tipped Applicator

Polystyrene handle



Actual Size

Description

This synthetic tipped applicator is made of medical-grade polyester fiber that is carded in our factory to maintain consistency. It is non-absorbent and contains no fatty acids or latent sugars, making it ideal for specimen transport, lab work, and diagnostic kits. It is bonded securely to a semi-flexible polystyrene handle by an aqueous based adhesive. Available sterile or non-sterile.

Dimensions

6" Overall Length (152.4mm)

Tip

■ 0.203" Dia. (5.16mm)

■ 0.687" Length (17.45mm)

Shaft

■ 0.100" Dia. (2.54mm)

■ 5.875" Length (149.23mm)

Standards

- CE marked for medical use
- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of compliance available with each shipment upon request

Qualities

- Medical grade quality
- Non-absorbent tip
- Contains no fatty acids or latent sugars

Applications

- Diagnostic specimen collection
- Patient care
- Cell collection for DNA testing
- Gynecological screening

Ref #	Packaging
806-PD <i>Non-sterile</i>	10/10/100 (10 boxes per case, 10 bags per box, 100 swabs per bag)
25-806 1PD <i>Sterile</i>	10/100/1 (10 boxes per case, 100 individually wrapped swabs per box, 1 per package)
25-806 2PD <i>Sterile</i>	10/100/2 (10 boxes per case, 100 individually wrapped swabs per box, 2 per package)



Quality Medical Products Since 1919

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

1-800-321-2313
puritanmedproducts.com