



Puritan
Quality since 1919

800-321-2313

PURITANMEDPRODUCTS.COM

DATA SHEET

Polyester Tipped Applicator - Dry Transport System

Puritan's polypropylene dry transport tube and plug-style cap is perfect for specimen collection and transportation to the lab or use in diagnostic kits. This synthetic tipped applicator is made of medical-grade polyester fiber that is carded on-site to maintain consistency. It is boned securely to a semi-flexible polystyrene handle by an aqueous based adhesive.

Dimensions

Overall Length

- Tube: 6.062" (153.975 mm)
- Swab: 6.0" (152.4 mm)

Swab Handle

- 0.1" Dia. (2.54 mm)
- 5.875" Length (149.225 mm)

Swab Tip

- 0.203" Dia. (5.156 mm)
- 0.687" Length (17.450 mm)

Qualities

- Medical grade quality
- Non-absorbent tip
- Contains no fatty acids or latent sugars
- Packaged in a reclosable dry transport tube for ease in transport to the lab
- Made in USA

Standards

- CE marked for medical use
- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of sterility and/or conformance available with each shipment upon request
- Made in USA

Applications

- Diagnostic specimen collection
- Forensic sampling
- Cell collection for DNA testing
- Gynecological screening
- Testing industrial surfaces



Ref#	Packaging	
25-806 1PD BT	5/100/1	5 boxes per case, 100 tubes per box, 1 swab per tube