

# Rayon Tipped Applicator

Polystyrene handle



- USP Grade Rayon Tip

- Semi-flexible polystyrene handle

Actual Size

## Description

This specialty swab is made with USP grade rayon fiber that is carded on-site to maintain consistency. This swab is ideal for diagnostic specimen collection, wound treatment and application of medicines. It is bonded securely to a semi-flexible polystyrene handle by an aqueous based adhesive.

## Dimensions

6" Overall Length (152.4 mm)

### Tip

■ 0.203" Dia. (5.156 mm)

■ 0.688" Length (17.475 mm)

### Handle

■ 0.100" Dia. (2.54 mm)

■ 5.875" Length (149.225 mm)

## Standards

- CE marked for medical use
- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of sterility and/or compliance available with each shipment upon request

## Qualities

- Medical grade quality
- Absorbent tip
- Soft smooth tip

## Applications

### Diagnostics / Forensics / Medical

- Specimen collection swab for testing & screening
- Cell collection for DNA testing
- General purpose collection swab

Ref #	Packaging
<b>25-806 1PR</b> <i>Sterile</i>	10/100/1 (10 boxes of 100 individually wrapped packages; 1 per package)
<b>25-806 2PR</b> <i>Sterile</i>	10/100/2 (10 boxes of 100 individually wrapped packages; 2 per package)



Quality Medical Products Since 1919

Puritan Medical Products Company LLC  
Guilford, Maine 04443-0149 USA

Tel 800-321-2313 (US and Canada Only)  
207-876-3311

Fax 800-323-4153  
207-876-3130

Web [www.puritanmedproducts.com](http://www.puritanmedproducts.com)

Email [sales@puritanmedproducts.com](mailto:sales@puritanmedproducts.com)

**1-800-321-2313**  
**[puritanmedproducts.com](http://puritanmedproducts.com)**