



Puritan
Quality since 1919

800-321-2313

PURITANMEDPRODUCTS.COM

DATA SHEET

Rayon Tipped Applicator - Dry Transport System

Polystyrene handle - sterile

Puritan's polypropylene dry transport tube and plug-style cap is perfect for specimen collection and transportation to the lab or use in diagnostic kits. This speciality applicator is made with USP grade rayon fiber that is carded on-site to maintain consistency. It is boned securely to a semi-flexible polystyrene handle by an aqueous based adhesive.

Dimensions

Overall Length

- Tube: 6.062" (153.975 mm)
- Swab: 6.0" (152.4 mm)

Swab Handle

- 0.1" Dia. (2.54 mm)
- 5.875" Length (149.225 mm)

Swab Tip

- 0.203" Dia. (5.156 mm)
- 0.687" Length (17.450 mm)

Qualities

- Medical grade quality
- Absorbent tip
- Soft and smooth tip
- Packaged in a reclosable dry transport tube for ease in transport to the lab
- Made in USA

Standards

- CE marked for medical use
- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of sterility and/or conformance available with each shipment upon request

Applications

- Specimen collection for testing and screening
- Cell collection for DNA testing
- General purpose collection swab



Ref#	Packaging	
25-806 1PR BT	5/100/1	5 boxes per case, 100 tubes per box, 1 swab per tube