



Puritan
Quality since 1919

800-321-2313

PURITANMEDPRODUCTS.COM

DATA SHEET

Calcium Alginate Tipped Applicator - Dry Transport System

Wood handle - sterile

Puritan's polypropylene dry transport tube and plug-style cap is perfect for specimen collection and transportation to the lab or use in diagnostic kits. This speciality applicator is tipped with calcium alginate fiber, an organic material derived from seaweed and carded in Puritan's facility to assure conformance with medical device specifications. It is boned securely to a rigid wood handle by an aqueous based adhesive.

Dimensions

Overall Length

- Tube: 6.062" (154 mm)
- Swab: 5.937" (150.8 mm)

Swab Handle

- 0.1" Dia. (2.54 mm)
- 5.812" Length (147.625 mm)

Swab Tip

- 0.219" Dia. (5.563 mm)
- 0.641" Length (16.281 mm)

Standards

- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of sterility and/or conformance available with each shipment upon request

Applications

- Specimen collection for testing and screening
- Wound care

Qualities

- Medical grade quality
- Calcium alginate tip contains no fatty acid, which may be inhibitory to many strains of bacteria
- Packaged in a reclosable dry transport tube for ease in



Ref#	Packaging	
25-806 1WA BT	5/100/1	5 boxes per case, 100 tubes per box, 1 swab per tube