



Puritan
Quality since 1919

800-321-2313

PURITANMEDPRODUCTS.COM

DATA SHEET

Cotton Tipped Applicator - Dry Transport System

Wood handle - sterile

Puritan's polypropylene dry transport tube and plug-style cap is ideal for specimen collection and transportation to the lab or for use in diagnostic kits. This absorbent cotton tipped applicator is produced with high-quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a sturdy wood handle by an aqueous based adhesive.

Dimensions

Overall Length

- Tube: 6.062" (154 mm)
- Swab: 5.937" (150.8 mm)

Swab Handle

- 0.1" Dia. (2.54 mm)
- 5.812" Length (147.625 mm)

Swab Tip

- 0.188" Dia. (4.775 mm)
- 0.625" Length (15.875 mm)

Standards

- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of sterility and/or conformance available with each shipment upon request

Applications

- Specimen collection for testing and screening
- Cell collection for DNA testing

Qualities

- High-quality USP grade cotton
- Packaged in a reclosable dry transport tube for ease in transport to the lab
- Made in USA



Ref#	Packaging	
25-806 1WC BT	5/100/1	5 boxes per case, 100 tubes per box, 1 swab per tube