

# Cap-Shure® Cotton Tipped Applicator

Wood handle



Aerated tip protector  
over cotton tip

Sturdy wood handle

\* Not actual size

## Description

This absorbent cotton tipped applicator is produced with High Quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a sturdy wood handle by an aqueous based adhesive. Aerated tip protector designed for convenient, reliable specimen, DNA or evidence collection. Comes as a unit with a reclosable, ventilated cap for air drying and transporting samples.

## Swab dimensions

6" Overall Length (152.4 mm)

### Tip

■ 0.188" Dia. (4.775 mm)

■ 0.625" Length (15.875 mm)

### Handle

■ 0.100" Dia. (2.54 mm)

■ 5.812" Length (147.625 mm)

## Standards

- CE marked for medical use
- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of sterility and/or compliance available with each shipment upon request

## Qualities

- Medical grade quality
- Tip suspended away from sides
- Can lay flat to dry

## Applications

### Diagnostic / Medical / Forensics

- Specimen collection
- DNA collection
- Evidence collection

Ref #	Packaging
<b>25-806 1WC EC</b> <i>Sterile</i>	10/50/1 (10 boxes per case; 50 individually wrapped packages per box; 1 per package)



Quality Medical Products Since 1919

Puritan Medical Products Company LLC  
Guilford, Maine 04443-0149 USA

**Tel** 800-321-2313 (US and Canada Only)  
207-876-3311

**Fax** 800-323-4153  
207-876-3130

**Web** [www.puritanmedproducts.com](http://www.puritanmedproducts.com)

**Email** [sales@puritanmedproducts.com](mailto:sales@puritanmedproducts.com)

**1-800-321-2313**  
**[puritanmedproducts.com](http://puritanmedproducts.com)**