



**Puritan**  
Quality since 1919

800-321-2313

PURITANMEDPRODUCTS.COM

## DATA SHEET

# Rayon Tipped Applicator - Dry Transport System

Wood handle - sterile

Puritan's polypropylene dry transport tube and plug-style cap is perfect for specimen collection and transportation to the lab or use in diagnostic kits. This speciality applicator is made with USP grade rayon fiber that is carded on-site to maintain consistency. It is boned securely to a rigid wood handle by an aqueous based adhesive.

## Dimensions

Overall Length

- Tube: 6.062" (153.975 mm)
- Swab: 5.937" (150.8 mm)

Swab Handle

- 0.1" Dia. (2.54 mm)
- 5.812" Length (147.625 mm)

Swab Tip

- 0.188" Dia. (4.775 mm)
- 0.625" Length (15.875 mm)

## Qualities

- Medical grade quality
- Absorbent tip
- Soft and smooth tip
- Packaged in a reclosable dry transport tube for ease in transport to the lab
- Made in USA

## Standards

- CE marked for medical use
- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of sterility and/or conformance available with each shipment upon request
- Made in USA

## Applications

- Specimen collection for testing and screening
- Cell collection for DNA testing
- General purpose collection swab



Ref#	Packaging	
25-806 1WR BT	5/100/1	5 boxes per case, 100 tubes per box, 1 swab per tube