

DNA-free Cotton Tipped Applicator

Wood shaft



- USP Grade Cotton Tip

- Northern Maine White Birch Shaft

Actual Size



Puritan®
Quality since 1919

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Description

This absorbent cotton tipped applicator is produced with High Quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a rigid white birch shaft by an aqueous based adhesive. Product is packed in an individually wrapped transport tube.

Tested by an independent lab to the most stringent standards (less than 25 pg of human DNA), certified to be DNA-free, DNase-free, and RNase-free.

Dimensions

6" Overall Length (152.4mm)

Tip

■ .188" Dia. (4.78mm)

■ .625" Length (15.88mm)

Shaft

■ .100" Dia. (2.54mm)

■ 5.812" Length (147.62mm)

Standards

- CE marked for medical use
- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certified to be DNA-free, DNase-free & RNase-free

Qualities

- Medical grade quality
- Absorbent tip
- Soft smooth tip

Applications

Forensics : DNA Testing

- Crime scene investigation
- Buccal cell collection
- Specimen collection

Ref #	Packaging
25-806 2WC TT FDNA <i>Sterile</i>	10/50/2 (10 boxes per case, 50 tubes individually wrapped per box, 2 swabs per tube)

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

1-800-321-2313
puritanmedproducts.com