


## Fab-Swab Cotton Tipped Applicator

### Wood handle



Not actual size

### Description

This DNA controlled absorbent cotton tipped applicator is produced with High Quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a rigid wood handle.

The convenient, reclosable, dry transport device is designed for reliable sample and evidence collection. The patented tube design entails a breathable filter label that covers holes to prohibit mold growth during transport. ID labels are included.

### Dimensions

5.095" Overall Length (129.413 mm)

#### Tip

- 0.12" Dia. tapered to 0.05"  
(3.048 mm tapered to 1.27 mm)
- 0.5" Length (12.7 mm)

#### Handle

- 0.083" Dia. (2.108 mm)
- 5.032" Length (127.813 mm)

### Standards

- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of compliance available with each shipment upon request

### Qualities

- Medical grade quality
- Produced in a DNA controlled environment
- Reclosable transport device
- Can be used dry, with IPA, or sterile water
- Includes evidence information label
- Tapered mini tip for reaching into small areas

### Applications

- DNA testing
- Forensic evidence collection

Ref #	Packaging	
<b>25-825 1WCS TT FABUSA</b> <i>Sterile</i>	10/50/1	(10 boxes per case; 50 individually wrapped tubes per box; 1 swab per per tube)



Quality Medical Products Since 1919

Puritan Medical Products Company LLC  
Guilford, Maine 04443-0149 USA

**Tel** 800-321-2313 (US and Canada Only)  
207-876-3311

**Fax** 800-323-4153  
207-876-3130

**Web** [www.puritanmedproducts.com](http://www.puritanmedproducts.com)

**Email** [sales@puritanmedproducts.com](mailto:sales@puritanmedproducts.com)

**1-800-321-2313**  
[puritanmedproducts.com](http://puritanmedproducts.com)