



Puritan
Quality since 1919

800-321-2313

PURITANMEDPRODUCTS.COM

DATA SHEET

ESK - Buffered Peptone Water

ESK® Sampling Kit - Pre-filled with 4ml or 10 ml Buffered Peptone Water

The Buffered Peptone Water is a pre-enriched nonselective buffered solution. It allows for the repair of damaged cells and facilitate the recovery of target bacteria in samples taken for sanitation monitoring in industrial, food and dairy plants.

Dimensions

Overall length of Kit: 4.06" (103.124 mm)

Swab Handle

- 0.100" Dia. (2.54 mm)
- 3.125" Length (79.375 mm)

Swab Tip

- 0.203" Dia. (5.156 mm)
- 0.687" Length (17.450 mm)

Qualities


- Easy edge-peel packaging
- Quick-turn leak-resistant cap
- High-quality Puritan polyester tipped swab
- Shelf-life - 18 months

Standards

- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of sterility and/or conformance available with each shipment upon request

Applications

- For environmental surface sampling
- Quality Control
- Detection of microorganisms



Ref#	Packaging
25-83004 PD BPW	50/1 (4ml - 50 individually wrapped kits with polyester swab)
25-83004 PDB BPW	6/50/1 (4ml - 6 boxes, 50 individually wrapped kits with polyester swab, 1 swab per package)
25-83010 PD BPW	50/1 (10ml - 50 individually wrapped kits with polyester swab)
25-83010 PDB BPW	6/50/1 (10ml - 6 boxes, 50 individually wrapped kits with polyester swab, 1 swab per package)