



Puritan
Quality since 1919

800-321-2313

PURITANMEDPRODUCTS.COM

DATA SHEET

ESK - Letheen Broth

ESK® Sampling Kit - Pre-filled with 4ml or 10 ml Letheen Broth with Tween 80

Letheen Broth is recommended for testing materials for preservatives where neutralization of disinfectants such as quaternary ammonium is important for antimicrobial activity.

Dimensions

Overall length of Kit: 4.06" (103.124 mm)

Swab Handle

- .100" Dia. (2.54 mm)
- 3.125" Length (79.375 mm)

Swab Tip

- 0.203" Dia. (5.156 mm)
- 0.687" Length (17.450 mm)

Qualities

- Easy edge-peel packaging
- Quick-turn leak-resistant cap
- High-quality Puritan polyester tipped swab
- Shelf-life - 18 months

Standards

- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of sterility and/or conformance available with each shipment upon request

Applications

- For environmental surface sampling
- Quality Control
- Detection of Microorganisms



Ref#	Packaging
25-83004 PD LB	50/1 (4ml - 50 individually wrapped kits with polyester swab)
25-83004 PDB LB	6/50/1 (4ml - 6 boxes, 50 individually wrapped kits with polyester swab, 1 swab per package)
25-83010 PD LB	50/1 (10ml - 50 individually wrapped kits with polyester swab)
25-83010 PDB LB	6/50/1 (10ml - 6 boxes, 50 individually wrapped kits with polyester swab, 1 swab per package)