



Puritan
Quality since 1919

800-321-2313

PURITANMEDPRODUCTS.COM

DATA SHEET

ESK - Neutralizing Buffer

ESK® Sampling Kit - Pre-filled with 4ml or 10 ml Neutralizing Buffer

An ideal choice for on-the-spot quality control checks. The Neutralizing Buffer is recommended for verification in your sanitation programs with equipment, according to the standard methods.

Dimensions

Overall length of Kit: 4.06" (103.124 mm)

Swab Handle

- 0.100" Dia. (2.54 mm)
- 3.125" Length (79.375 mm)

Swab Tip

- 0.203" Dia. (5.156 mm)
- 0.687" Length (17.450 mm)

Qualities


- Easy edge-peel packaging
- Quick-turn leak-resistant cap
- High-quality Puritan polyester tipped swab
- Shelf-life - 18 months

Standards

- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of sterility and/or conformance available with each shipment upon request

Applications

- For environmental surface sampling
- Quality Control
- Detection of Microorganisms



Ref#	Packaging	
25-83004 PD NB	50/1	(4ml - 50 individually wrapped kits with polyester swab)
25-83004 PDB NB	6/50/1	(4ml - 6 boxes, 50 individually wrapped kits with polyester swab, 1 swab per package)
25-83010 PD NB	50/1	(10ml - 50 individually wrapped kits with polyester swab)
25-83010 PDB NB	6/50/1	(10ml - 6 boxes, 50 individually wrapped kits with polyester swab, 1 swab per package)