



Puritan
Quality since 1919

800-321-2313
PURITANMEDPRODUCTS.COM

DATA SHEET

Environmental Sampling Swab

EnviroMax® Round dry foam tip - Sterile or Non-Sterile

Ready to use environmental foam sampling swabs, provided in wet and dry format for accuracy and consistency of sample collection.

Dimensions

Overall length: 6.02" (152.908 mm)

Handle

- 0.125" Dia. (3.175 mm)
- 5.60" Length (142.240 mm)

Tip

- 0.612" Dia. (15.545 mm)
- 2.220" Length (56.388 mm)

Qualities

- Large Tip for maximized collection and elution
- Rigid foam paddle tip
- High-quality Puritan Foam tipped swab
- Shelf-life - 12 months

Standards

- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of sterility and/or conformance available with each shipment upon request

Applications

- For environmental surface sampling
- Perfect for Spot Checks
- Quality Control
- Detection of Microorganisms



Ref#	Packaging	
88010-PF RND	25/1	(25 individually wrapped tubes with foam swab)
88010-PFB RND	10/25/1	(10 boxes, 25 individually wrapped kits with foam swab, 1 swab per package)
25-88010 PF RND	25/1	(25 individually wrapped tubes with foam swab)
25-88010 PFB RND	10/25/1	(10 boxes, 25 individually wrapped tubes with foam swab, 1 swab per package)



EnviroMax®