

Knitted Polyester Applicators

Static dissipative polypropylene handle



- Flexible paddle tip

- Solid static-dissipative polypropylene handle

Actual Size

Description

This knitted polyester tipped applicator is non-abrasive and lint-free. The tip can easily stand up to IPA, acetone, and other solvents commonly used in clean rooms and other critical environments. It is ultra-sonically welded to a solid static dissipative handle.

Dimensions

6.018" Overall Length (152.857 mm)

Tip

■ 0.661" Length (16.789 mm)

■ 0.268" Dia. (6.807 mm)

Handle

■ 5.922" Length (150.418 mm)

■ 0.12" Dia. (3.048 mm)

Standards

■ ISO 9001:2008

■ FDA GMP's

Qualities

- Lint and residue free
- Soft and non-abrasive
- Can be used with IPA and DI water, or dry
- Excellent wipe-dry characteristics
- Chemical resistant

Applications

Critical Environment/Cleanroom/Industrial

- Removing adhesives or other fluids
- Applying solutions or solvents
- General cleaning
- Used in electronic assembly and medical device manufacturing

Ref #	Packaging
3676 ESD	20/50 (20 resealable poly bags per case, 50 swabs per bag)



Specialty Applicators

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

1-800-321-2313
puritanmedproducts.com