


HydraFlock® Collection Device

Elongated tip, polystyrene handle



Description

Engineered to provide maximum absorption and elution of specimens. The unique, proprietary microstructure of HydraFlock's fibers enhances rapid absorption and release of biological specimens.

Patented: US 8334134, Canada 2,744,282

Patent pending: Europe 10251388.4

Dimensions

6" Overall Length (152.4 mm)

Tip

■ 0.219" Dia. (5.562 mm)

■ 0.685" Length (17.399 mm)

Handle

■ 0.1" (2.54 mm) tapered to 0.067" (1.702 mm)

■ 5.96" Length (151.384 mm)

Standards

- CE marked for medical use
- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of sterility and/or conformance available with each shipment upon request

Qualities

- Medical grade quality
- Highly absorbent tip
- Elongated, soft smooth tip

Applications

- For Buccal cell collection
- Rapid diagnostic testing
- Specimen collection
- Forensic evidence collection
- DNA collection

Ref #	Packaging
3806-H <i>Non-sterile</i>	10/100 (10 bags per case, 100 per bag)
25-3806-H <i>Sterile</i>	10/50/1 (10 boxes per case, 50 individually wrapped packages, 1 swab per package)



HydraFlock®

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

1-800-321-2313
puritanmedproducts.com