



Puritan
Quality since 1919

800-321-2313

PURITANMEDPRODUCTS.COM

DATA SHEET

Self-Saturating Applicator

Popule® Self-saturating foam tip, handle filled with 91% Isopropyl alcohol and 9% distilled water.

This pre-filled popule is made with a rectangular 100PPI medical grade foam tip. It is bonded to a clear polypropylene tube. Snap the handle to release the liquid into tip - neat and clean for your convenience!

Dimensions

Overall length: 4.125" (104.775 mm)

Handle

- 0.200" Dia. (5.08 mm)
- 3.99" Length (101.52 mm)

Foam Tip

- 0.562" Dia. (14.275 mm)
- 1.062" Length (26.975 mm)

Qualities


- Lint and residue free
- Non-glass handle
- Combines solvent and applicator in one
- No cross contamination
- PurSwab brand Popule
- Made in U.S.A

Standards

- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of sterility and/or conformance available with each shipment upon request

Applications

- Ideal for cleaning keyboards and other electronic equipment
- General removal of residues
- Used in critical environment applications



Ref#	Packaging
4545	20/50 (20 bags, 50 swabs per bag)



Puritan®