

Rayon Tipped Applicator

Polystyrene handle



- USP Grade Rayon Tip

- Semi-flexible polystyrene handle

Actual Size

Description

This specialty swab is made with USP grade rayon fiber that is carded on-site to maintain consistency. This swab is ideal for diagnostic specimen collection, wound treatment and application of medicines. It is bonded securely to a semi-flexible polystyrene handle by an aqueous based adhesive.

Dimensions

6" Overall Length (152.4 mm)

Tip

■ 0.203" Dia. (5.156 mm)

■ 0.687" Length (17.45 mm)

Handle

■ 0.100" Dia. (2.54 mm)

■ 5.875" Length (149.225 mm)

Standards

■ CE marked for medical use

■ ISO 9001:2008

■ ISO 13485:2003

■ FDA GMP's

Qualities

■ Medical grade quality

■ Semi-flexible handle

■ Soft smooth tip

Applications

Medical / Health-Beauty / Industrial / Veterinary

■ General purpose applicator

■ Specimen collection

■ Application of topical dressings & first aid

■ Industrial general cleaning

■ Cleaning sensors and components

■ Removing residues

Ref #	Packaging
806-PR <i>Non-sterile</i>	10/10/100 (10 boxes per case, 10 bags per box, 100 swabs per bag)



Quality Medical Products Since 1919

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

1-800-321-2313
puritanmedproducts.com