


Cotton Tipped Applicator

Wood handle, No glue



Description

This absorbent cotton tipped applicator is produced with High Quality USP grade cotton which is carded and coiled onsite to maintain consistency. It has a low-lint tip that is ideal for cleaning lint sensitive components.

Dimensions

3.032" Overall Length (77.013 mm)

Tip

■ 0.156" Dia. (3.962 mm)

■ 0.625" Length (15.875 mm)

Handle

■ 0.083" Dia. (2.108 mm)

■ 2.906" Length (73.812 mm)

Standards

■ ISO 9001:2008

■ ISO 13485:2003

■ FDA GMP's

■ Certificate of compliance available with each shipment upon request

Qualities

■ Small cotton tip, low- lint

■ Sturdy wood handle

■ Made in USA

Applications

Critical Environment/Industrial/Veterinary

■ Used in general cleaning of smaller areas

■ Ideal for cleaning lint sensitive items

Ref #	Packaging
8380-WCS <i>Non-sterile</i>	10/10/100 (10 boxes per case, 10 bags, per box, 100 swabs per bag)



Puritan®

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

1-800-321-2313
puritanmedproducts.com