

Cotton Tipped Applicator

Semi-flexible polystyrene handle

Description

This absorbent double cotton tipped applicator is produced with High Quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to an semi-flexible plastic handle by an aqueous based adhesive.

Dimensions

6.125" Overall Length (155.575 mm)

Tip

■ 0.203" Dia. (5.156 mm)

■ 0.687" Length (17.45 mm)

Handle

■ 0.100" Dia. (2.54 mm)

■ 5.875" Length (149.225 mm)

Standards

■ ISO 9001:2008

■ ISO 13485:2003

■ FDA GMP's

■ Certificate of compliance available with each shipment upon request

Qualities

■ Standard cotton-tip/ Double

■ Semi-flexible hollow polystyrene handle

■ Made in USA

Applications


Critical Environment/Health&Beauty/Medical

■ All purpose general swab

■ For use in rework and repair

■ Removing excess solder and fluxes from circuit boards

■ Cleaning metal finishes



Ref #	Packaging
861-PC DBL <i>Non-sterile</i>	10/5/100 (10 boxes per case, 5 bags per box, 100 swabs per bag)



Puritan®

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

1-800-321-2313
puritanmedproducts.com