


Cotton Tipped Applicator

Double cotton tip, Wood handle



Description

This absorbent double cotton tipped applicator is produced with High Quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a wood handle by an aqueous based adhesive.

Dimensions

6.062" Overall Length (153.975 mm)

Tip

■ 0.188" Dia. (4.775 mm)

■ 0.625" Length (15.875 mm)

Handle

■ 0.100" Dia. (2.54 mm)

■ 5.812" Length (147.625 mm)

Standards

■ ISO 9001:2008

■ ISO 13485:2003

■ FDA GMP's

■ Certificate of compliance available with each shipment upon request

Qualities

■ Standard cotton-tip/ Double

■ Sturdy wood handle

■ Made in USA

Applications

Critical Environment/Health&Beauty/Medical

■ All purpose general swab

■ For use in rework and repair

■ Removing excess solder and fluxes from circuit boards

■ Cleaning metal finishes

Ref #	Packaging
861-WC DBL <i>Non-sterile</i>	10/5/100 (10 boxes per case, 5 bags per box, 100 swabs per bag)



Puritan®

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

1-800-321-2313
puritanmedproducts.com