

Cotton Tipped Applicator

Polystyrene handle



- USP Grade Cotton, small tip

- Semi-flexible polystyrene handle

Actual Size

Description

This cotton tipped applicator is produced with High Quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a semi-flexible polystyrene handle by an aqueous based adhesive. Small cotton tip ideal for hard-to-reach or small areas.

Dimensions

■ 6" Overall Length (152.4 mm)

Tip

■ 0.172" Dia. (4.369 mm)

■ 0.687" Length (17.45 mm)

Handle

■ 0.1" Dia. (2.54 mm)

■ 5.875" Length (149.225 mm)

Standards

- FDA GMP's
- ISO 9001:2008
- ISO 13485:2003

Qualities

- Small tip for hard-to-reach or small areas
- Semi-flexible polystyrene handle
- Features minimal cotton for small or hard-to-reach areas

Applications

Critical Environment / Medical

- General cleaning swab
- Removing oil or dirt from machine parts

Ref #	Packaging
868-PCS <i>Non-sterile</i>	10/10/100 (10 boxes per case; 10 bags per box; 100 swabs per bag)



Quality Medical Products Since 1919

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

1-800-321-2313
puritanmedproducts.com