



Puritan
Quality since 1919

800-321-2313
PURITANMEDPRODUCTS.COM

DATA SHEET

Cotton Tipped Applicator

Lint-less, double pointed tip, fine paper handle.

This absorbent cotton double-tipped applicator is produced with pharmaceutical grade cotton. It is bonded securely to a paper shaft by an aqueous based adhesive.

Dimensions

3.11" Overall Length (78.994 mm)

Handle

- 0.060" Dia. (1.524 mm)
- 2.874" Length (73.0 mm)

Tip

- 0.118" Dia. (2.997 mm)
- 0.433" Length (10.998 mm)

Qualities

- Lint-less and residue free
- Easy disposal
- Highly absorbent cotton does not swell in solvents
- Made in USA

Standards

- CE marked for medical use
- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of sterility and/or conformance available with each shipment upon request

Applications

- Cleans sensitive components
- Cleaning disk drives, optics, and other critical components



Ref#	Packaging
871-PC DBL (Non-Sterile)	10/25 (100 leaves per pouch, 25 swabs per leaf)



Puritan®