

Cotton Tipped Applicator

Anti-static Polypropylene handle



Description

This absorbent double cotton tipped applicator is produced with High Quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to an anti-static polypropylene handle by an aqueous based adhesive.

Dimensions

6" Overall Length (152.4 mm)

Tip

■ 0.203" Dia. (5.156 mm)

■ 0.688" Length (17.475 mm)

Handle

■ 0.100" Dia. (2.54 mm)

■ 5.875" Length (149.225 mm)

Standards

■ ISO 9001:2008

■ ISO 13485:2003

■ FDA GMP's

■ Certificate of compliance available with each shipment upon request

Qualities

■ Lint free

■ Anti-static polypropylene handle

■ Regular cotton tip

■ Made in USA

Applications

Critical Environment/Industrial

■ Cleaning and removing residues

Ref #	Packaging
876-PPC <i>Non-sterile</i>	10/10/100 (10 boxes per case, 10 bags per box, 100 swabs per bag)



Puritan®

Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

1-800-321-2313
puritanmedproducts.com