



DATA SHEET

Cotton Tipped Applicator

Lint-free cotton tip, wood handle.

This absorbent cotton tipped applicator is produced with high quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a naturally biodegradable wooden handle with a hot-melt glue adhesive.

Dimensions

5.937" Overall Length (150.799 mm)

Handle

- 0.100" Dia. (2.54 mm)
- 5.812" Length (147.625 mm)

Tip

- 0.188" Dia. (4.775 mm)
- 0.625" Length (15.875 mm)

Standards

- CE marked for medical use
- ISO 9001:2008
- ISO 13485:2003
- FDA GMP's
- Certificate of sterility and/or conformance available with each shipment upon request

Qualities

- Medical grade quality
- Highly absorbent tip
- Lint-free
- Soft smooth tip

Applications

- General cleaning swab
- Cleaning cameras
- Sensitive components requiring low lint

Ref#	Packaging
876-WC (Non-sterile)	10/10/100 (10 boxes per case, 10 bags per box, 100 per bag)



Puritan®