

896-PC PurSwab Cotton Tipped Applicator



Features and Benefits

Tip:	Shaft:	Uses:	Packaging:	Qualities:	Extraction Analysis*
• Single-end cotton • Smooth • Cylindrical • 4mm diameter • 12mm length	• Paper • White • 146mm length • 2.46mm diameter	• Health and beauty • Cleaning disk drives, optics, PC board assembly and other critical industrial components • Cleans sensitive components • Veterinary	• 1,000 per case (50 pcs packed in a heat sealed poly bag, 20 bags of 50 in a case, 20/50) • No silicone • Anti-static packaging	• Easy disposal • Lint-less and residue free • Highly absorbent cotton does not swell in solvents • Low in NVR • Burnable • Environmental friendly	• Chloride (Cl) 7.33 µg/tip • Sulphate 7.54 µg/tip • Flouride 0.68 µg/tip • Sodium 10.85µg/tip • NVR (Di) 0.25 mg/tip • NVR (IPA) 0.11 mg/tip • NVR (Acetone) 0.13 mg/tip • NVR (IPA:Hexane) 0.17 µg/tip • Silicone oil N/D mg/tip * Based on a typical sample batch



Puritan Medical Products Company LLC
Guilford, Maine 04443-0149 USA

Tel 800-321-2313 (US and Canada Only)
207-876-3311

Fax 800-323-4153
207-876-3130

Web www.puritanmedproducts.com

Email sales@puritanmedproducts.com

1-800-321-2313
puritanmedproducts.com