



- Separates into two pieces for easy carrying
- Rotation of pin either locks or allows ramp to be separated into two pieces (shown below)
- Folds in half, carries like a suitcase
- Patented hinge eliminates pinch point
- Closure strap locks ramp panels together
- Lightweight, easy-to-handle & set up
- Durable welded construction
- Safety guide indicates if ramp is on a safe slope
- Full platform provides excellent stability
- Anti-slip, high traction surface
- Accommodates wheelchairs and scooters with various wheel configurations
- Safety DVD and steel security pins included
- This ramp has been tested at a 3X safety factor



**Custom widths are available. Give your specifications to your local PVI retailer!**

Item Number	Length	Width	Folded Size	Weight	Capacity*
WCR530	5'	30"	32¼" x 16" x 9"	30 lbs.	800 lbs.
WCR630	6'	30"	38¼" x 16" x 9"	35 lbs.	800 lbs.
WCR730	7'	30"	44¼" x 16" x 9"	39 lbs.	800 lbs.
WCR830	8'	30"	50¾" x 16" x 9½"	52 lbs.	800 lbs.
WC1030	10'	30"	62¾" x 16" x 9½"	63 lbs.	800 lbs.
WC1230	12'	30"	74¾" x 16" x 9½"	73 lbs.	600 lbs.

**LIMITED LIFETIME WARRANTY**

\* Maximum ramp capacity: 400 lbs. (one axle) – 800 lbs. (two axles). Except on WC1230. Maximum ramp capacity: 300 lbs. (one axle) – 600 lbs (two axles).

