

# QuickVue® Influenza Nasal Swab Transport Tubes

25 foam-tipped nasal swabs inside plastic transport tubes.

For use with QuickVue Influenza (Cat. #00317) and QuickVue Influenza A+B (Cat. #20183) tests.  
Cat. #20171



Specimens should be tested as soon as possible after collection. For appropriate storage conditions prior to testing, please refer to the Specimen Collection and Handling section in the package insert.

**Note:** If dual testing for both influenza and RSV using the QuickVue Influenza A+B and RSV test kits, nasal swabs should not be used. Nylon flocked nasopharyngeal swabs made by Copan Diagnostics are recommended for use. Nasopharyngeal (NP) swabs may be ordered using Quidel Catalog #20226 (NP swabs only). These NP swabs are also available in dry transport tubes and may be ordered through Diagnostic Hybrids Inc., or various distributors, using Copan Item #551C.



To order: Call your Quidel Account Manager at 1.800.874.1517. Outside the U.S. call: 858.552.1100.

QuickVue Swab Transport Tube Set  
Catalog #20171



For more information on QuickVue products, visit our website or call your Quidel Account Manager at 1.800.874.1517 (U.S. only), or 858.552.1100 (outside U.S.).

Quidel Corporation • Worldwide Headquarters • 10165 McKellar Court • San Diego, CA 92121 USA • [quidel.com](http://quidel.com)

Distributor Information: