

SPECIFICATION

Table 2

ITEM	UNIT	888WA
Overall Length ¹	Inches	38"
Overall Width	Inches	22"
Overall Height	Inches	40"
Total Weight (w/ 26 Ah Batteries)	Lbs	160
Battery Weight	Lbs*pcs	17 lbs*2
Battery Capacity	12V_Ah x 2pcs	26Ah
Safe Climbing Angle	°(%)	8(14)
Suggested User Weight (on level surface)	Lbs	Max. 300
Maximum Speed ³	mph	5
Range ⁴ (Per charge w/ standard 26 Ah battery)	Miles	8
Charger Type	On / Off Board	Off Board 24V, 5A

All specifications are subject to change without prior notice. SHOPRIDER® / PIHSIANG® reserves the

rights of any changes on the unit.

1. Includes the anti-tip wheel.

2. The actual driving range varies with the factors shown below.

(1) the weight of occupant

(2) ground surface

(3) battery capacity and conditions

(4) type of charger

(5) ambient temperature

(6) Driving methods

(7) If the battery and mechanical moving parts fully broken in.