

Scootie™

(TE-787NA)

Scootie™ was the first to set the standard in compact scooters that no other brand can match. Light-weight, easy to transport, articulating front wheels...are just some of the features that the **Scootie™** has. The **Scootie™** is only 82 lbs. with batteries. The patented articulating front wheels allow four-wheel stability with three-wheel maneuverability. The **Scootie™** is length adjustable and can fold down for easy transport.



Available Colors



Red



Blue

- 250 lb. weight capacity
- Adjustable flip-up arms
- Length adjustable
- Removable battery pack
- Puncture proof tires
- Connector-less battery pack

Overall Dimensions (L' x W x H)	in	37-40 x 17 x 34
Number, Size of Tires (front; rear)	Front (rear)	2, 7" (2, 8")
Suggested User Weight (on a level road)	lb	Max. 250
Weight of Heaviest Part (w/o baskets and seat)	lb	25
Total Weight (w/ batteries, basket, seat)	lb	82
Battery Pack Weight	lb	24
Battery Module Capacity	_V_Ah x Pcs	12V12Ah x 2
Charger		Off Board (2A)
Output Power of Motor	hp	.75
Maximum Speed ²	mph	3.5
Max. Climbing Angle/Safe Climbing Angle (based on user weight of 90kg/200lb.)	% (°)	14 (8)
Range (per charge) (After the battery and mechanical moving parts fully break in)	mile	10
Turning Radius	in	38
Ground Clearance	in	2
Handle Bar	Type	T-Bar

All specifications are subject to change without prior notice. Shoprider Mobility Products, Inc. reserves the rights of any changes to the unit.

1) Include the anti-tip wheel or the rear castor.

2) Driver weight may exceed weight of the unit; speed must be reduced when turning.

3) The actual driving range varies with the factors shown as below:

a) The weight of occupant

b) Ground surface

c) Battery capacity and conditions

d) Type of charger

e) Ambient temperature

f) The way of driving

g) If the battery and mechanical

moving parts fully break in

h) Etc.

SHOPRIDER®

Toll-Free: 800.743.0772

www.shoprider.com