



# Minimize the ***OUCH!*** Terumo® – The Sharper Needle



 **TERUMO®**

# Tests confirm

## What Veterinarians who use Terumo needles have known all along

Terumo needles are known worldwide for extreme sharpness. Our double-bevel design is precision-ground and honed to exacting quality standards. This attention to detail enables Terumo needles to penetrate skin with less force, which can mean greater comfort for patients and less stress for pet owners.


## Don't take our word for it...see it for yourself

A recent study\* confirms what veterinarians who use Terumo needles have known all along – Terumo needles are sharper than the biggest-selling veterinary brand. And as you and your patients know, sharper is better.

## Less penetration force means Sharper needles

	Penetration Force (gms) for Selected Needle Gauges			
	20G	22G	25G	30G
Terumo	11.0	10.2	8.8	6.7
Vet Market Leader	14.4	12.5	9.4	8.9

**On average,  
Terumo is  
18% sharper  
than the Vet  
market leader\***



\*Data on file.

**So try Sur-Vet for yourself!**

**Ask your representative** for samples of **SUR-VET** needles, and discover their superior sharpness for yourself.

**From the makers of SurFlo® I.V. Catheters**



TERUMO MEDICAL CORPORATION • 2101 Cottontail Lane • Somerset, New Jersey, 08873 • 800.283.7866 • [www.terumomedical.com](http://www.terumomedical.com)