



FLASH YOUR CONFIDENCE

With SurFlash[®] safety catheters—
clinically proven to increase first-stick success by up to 25%!¹



Because every stick counts.

Introducing SurFlash® safety catheters— a better experience for you and your patients

- The only line of safety IV catheters with a double-flash feature that verifies both needle **and** catheter placement
 - Reduced risk of back vein wall puncture
 - Clinically proven to increase first-stick success
 - Improved patient experience and satisfaction
 - Cost savings/fewer wasted catheters
- Passive safety
 - Shown to reduce needlestick injuries by 50%–70% over nonpassive safety²
- Full needle coverage; not just the needle tip
 - Reduced risk of blood contact

A proprietary safety valve puts the *Plus* in SurFlash® Plus

- Minimizes the risk of blood exposure between needle removal and IV attachment for a cleaner dressing site and IV hub area
- Reduces physical and visual exposure to blood

With *all* these advantages, SurFlash® and SurFlash® Plus are a dramatic upgrade over anything that has come before. Which leaves just one question: *Why would you use anything else?*



Flashes of insight

"We are very interested in a product that is clinically proven to increase first stick percentage."
—94% of nurses in a recent survey

Based on a survey of 51 nurses:

- 98% of nurses think a successful first stick is very important for patients' satisfaction
 - 96% say it's important for their *own* satisfaction
 - 98% believe the need for a second attempt at IV insertion increases the patient's anxiety
- 94% of nurses are very interested in a product that is clinically proven to increase first-stick percentage by up to 25%
- 94% also prefer an IV catheter that provides a visual indicator when the catheter is in the vessel, not just when the needle tip enters the vessel
- 90% would prefer to minimize the risk of blood exposure during IV insertions without the need to apply venous compression

**The bottom line is that nurses know what they want.
And SurFlash® has it all.**



Because every stick counts.

A flashy profile



Together, SurFlash® and SurFlash® Plus deliver the features and benefits nurses value most...

- Unique double-flash that verifies both needle **and** catheter placement
- Increased first-stick success
- Passive safety
- Reduced exposure to blood and needlestick injuries

...and more

- Pressure rating (up to 325 PSI)
- Bevel-up indicator
- Completely covered cannula
- Ultrasharp needle with medical-grade lubricant for smoother venipuncture
- Ergonomic shape for superior comfort and deployment



Because every stick counts.

SurFlash®

50/box, 4 boxes/case (200)		
Product Code	Needle Gauge	Needle Length
SR*FNP1832	18G	1 1/4"
SR*FNP2025	20G	1"
SR*FNP2032	20G	1 1/4"
SR*FNP2225	22G	1"
SR*FNP2232	22G	1 1/4"
SR*FNP2419	24G	3/4"

SurFlash® Plus

50/box, 4 boxes/case (200)		
Product Code	Needle Gauge	Needle Length
SR*FVP1832	18G	1 1/4"
SR*FVP2025	20G	1"
SR*FVP2032	20G	1 1/4"
SR*FVP2225	22G	1"
SR*FVP2232	22G	1 1/4"
SR*FVP2419	24G	3/4"

So go ahead. Put your hands on a SurFlash® safety catheter. And flash *your* confidence. For ordering information, visit terumotmp.com or call 1-800-888-3786.

References: 1. Data on file. 2. Tosini W, Ciotti C, Goyer F, et al. Needlestick injury rates according to different types of safety-engineered devices: results of a French multicenter study. *Infect Control Hosp Epidemiol.* 2010;31(4):402-407.

©, TERUMO, and SurFlash are trademarks owned by Terumo Corporation, Tokyo, Japan, and are registered with the U.S. Patent and Trademark Office.
©2012 Terumo Medical Corporation 11/12
All rights reserved. Accession 12072012-2520



TERUMO MEDICAL PRODUCTS