EasyTravel Scooter

Classic



Advanced Technologies That Move You

Combining the stable ride of 4 wheels with excellent maneuverability, this extremely lightweight scooter can be taken apart and folded in a few seconds.















CAR

IN

AIRPI ANE

STAIRS

FLEVATOR

EasyTravel Classic®

The total weight of the Classic is only 28kg (62 lbs) and separated into 3 parts, the heaviest part weighs a mere 13kg (29 lbs). With no connectors, latches or mechanisms to open when folding or taking apart, the Classic is simply the ideal scooter for all your running around.

Technical Specifications:		
Carrying Capacity	115 kg	250 lbs
Dimensions:		
Length	108 cm	43 in
Width	57 cm	22 in
Height	86 cm	34 in
Folded (vertical)	80x57x49 cm	31x22x19 in
Weight:		
Total Weight	28 kg	62 lbs
Front Column	9 kg	20 lbs
Rear Frame	13 kg	29 lbs
Battery:		
7 A/h	6 kg	13 lbs
12 A/h	8 kg	18 lbs
Range:		
7 A/h battery	up to 8 km	up to 5 miles
12 A/h battery	up to 13 km	up to 8 miles

Weights & measurements are approximate

- Quick and powerful, yet safe and stable.
- Folds easily and quickly. No need to bend or lift heavy parts.
- Ergonomically designed for comfort and functionality.
- Carry the Classic in your car, on a train or check into a plane, it travels anywhere with you.
- Detach the battery pack and carry it inside for recharging.
- Great for indoors and outdoors. Non-marking and flat free tires.





When folded, you can pull the Slassic like a suitcase



The Classic stored in the trunk of a car



The heaviest part weighs a mere 13 kg (29 lbs)









FOLD

