

EasyTravel Scooter

Elite[®]



Advanced Technologies That Move You

Simple, yet ingenious, the lightweight Elite[®] scooter folds in seconds without effort, bending or tools.



CAR



TRAIN



SHIP



AIRPLANE



STAIRS



ELEVATOR

EasyTravel Elite®

With no connectors, latches or mechanisms to open when folding or taking apart, the Elite can be easily placed in the back of a car, at your side on a train or a bus, or checked into a plane. **The Elite travels anywhere!**

Technical Specifications:

Carrying Capacity	115 kg	250 lbs
-------------------	--------	---------

Dimensions:

Length	108 cm	43 in
--------	--------	-------

Width	57 cm	22 in
-------	-------	-------

Height	93 cm	37 in
--------	-------	-------

Folded (vertical)	80x57x53 cm	31x22x21 in
-------------------	-------------	-------------

Weight:

Total Weight	29 kg	64 lbs
--------------	-------	--------

Front Column	11 kg	24 lbs
--------------	-------	--------

Rear Frame	12 kg	27 lbs
------------	-------	--------

Battery:

7 A/h	6 kg	13 lbs
-------	------	--------

12 A/h	8 kg	18 lbs
--------	------	--------


Range:

7 A/h battery	up to 8 km	up to 5 miles
---------------	------------	---------------

12 A/h battery	up to 13 km	up to 8 miles
----------------	-------------	---------------

Weights & measurements are approximate

- The total weight of the Elite is only 29Kg (64 Lbs) and separated into 3 parts, the heaviest part weighs a mere 12Kg (27 Lbs).
- Quick and powerful, yet safe and stable.
- Folds easily and quickly. No need to bend or lift heavy parts.
- Ideal for running errands or for just wandering around.
- Ergonomically designed for comfort and functionality.
- Detach the battery pack and carry it inside for recharging.
- Great for indoors and outdoors. Non-marking and flat free tires.
- Patented Hub Motor Drive Wheel.



When folded, you can pull the Elite like a suitcase



The Elite stored in the trunk of a car

Studio Orbit Jeger



SPLIT

The heaviest part weighs a mere 12Kg (27 Lbs)



FOLD