EasyTravel Mid Range Scooter Titan®



Advanced Technologies That Move You

The Titan® provides a unique combination of performance, comfort and portability. Experience its speed, comfort ride suspension and ease of transport.















EasyTravel Mid Range Titan®

Newly designed, ergonomic, portable & foldable, indoor & outdoor electric scooter

- Allows you to keep an independent, spontaneous and active life.
- Foldable and detachable in seconds without effort or bending!
- Compact, highly maneuverable, allowing it to be transferred easily into a car, a bus, a train or a plane.
- Ergonomic location of battery pack.

| Dimensions: | | |
|--------------------|--------|--------------|
| Length | 139 cm | (55 in) |
| Width | 61 cm | (24 in) |
| Height | 97 cm | (38 in) |
| Total Weight | 64 kg | (141 lb) |
| 2 x Battery 26 A/h | 9.5 kg | (21 lb) each |
| Front Frame | 22 kg | (48 lb) |
| Rear Frame | 23 kg | (51 lb) |
| Carrying Capacity | 136 kg | (300 lb) |
| Range up to | 30 km | (19 miles) |











