

Compact, All-Steel Construction,

Foot Operated Lid

- 24 Quart Capacity
- Flame Retardant
- Baked Epoxy Paint Finish
- All-Steel Construction
- Heavy-Duty Foot Pedal
- Self Closing Lid
- Easy To Clean
- Resists Scratching and Stains
- 1 Year Limited Warranty



The UMF Medical 1473R 24Qt. (22.7L) Waste Receptacle is made of heavy-duty steel and includes a baked epoxy paint finish for added durability. The foot operated lid is easily used through a heavy-duty foot pedal. It also includes a bag securing mechanism for quick and easy clean up, as well as free biohazard and infectious waste adhesive labels for OSHA labeling requirements. These waste receptacles are flame retardant.

The UMF Medical 1473R Waste Receptacle comes with a 1 year limited warranty.

Standard Features

- Flame retardant
- Unitized all-steel construction
- Baked epoxy paint finish for added durability against scratching and stains
- Full width hinge
- Heavy-duty foot pedal and linkage
- Self-closing lid
- Exposed edges hemmed for safety
- Bag securing mechanism
- Free biohazard and infectious waste adhesive labels for OSHA labeling requirements
- 1 year limited warranty

Specifications

- **Height:** 16" (40.6cm)
- **Width:** 11.75" (29.8cm)
- **Depth:** 13" (33.0cm)

Contact:

UMF Medical, 1316 Eisenhower Boulevard, Johnstown, PA 15904-3393
Phone: 814-266-8726, Toll Free: 800-638-5322, Fax: 814-266-1870 www.umfmedical.com