

5570

5500 Series Treatment Tables



Treatment Table with Adjustable Backrest



*Shown with optional accessory 015 – Foot section and pad.

Standard Features

- 500lb (226.8kg) weight capacity
- Ultra comfort seamless top with antimicrobial guard
- Manual, Adjustable Backrest
- Available in all upholstery color options
- All-welded steel construction
- Leveling glides

Optional Accessories

- **015** – Foot Section and pad
- **800** – Paper cutter band
- **804** – Paper rod holder
- **826** – Wood laminate storage shelf
- **5570F**– Fire retardant foam

Table Specifications

- **Height:** 31.5" (80cm)
- **Length:** 70"- 81" with leg section (177.8cm - 205.74cm)
- **Width:** 26.5" (67.3cm)

Colors

