

Phlebotomy Chair with Foot Operated Pump

*PATIENT COMFORT
Built in pillow headrest*

*ACCOMMODATING
Hinged armrests*

*INFECTION CONTROL
Seamless upholstery
with antimicrobial
guard*

*ACCESSIBILITY
360 degree rotation*

*PATIENT SAFETY
Locking swivel*

*EFFICIENCY
Interlinked
backrest and leg
section*

*INFECTION CONTROL
Bacteria fighting base
paint option*

Color Chart and Part Number:

8675RI River Rock	8675SN Sand Grey	8675MO Mocha	8675AD Adobe
8675ST Steel Blue	8675MI Midnight Blue	8675SG Sage	8675ON Onyx

Standard Features

- Sturdy steel frame
- Ultra comfort seamless top with antimicrobial guard
- Available in all UMF Medical upholstery options
- 360 degree rotation with swivel lock
- 10 degree elevation from horizontal position
- Manual foot pedal height adjustment
- Interlinked backrest and leg section for reclining patient to a horizontal position
- Hinged arms
- 250 pound (113.4kg) patient weight capacity
- 3 year limited warranty

Specifications

- **Length:** 67" (170.2cm) fully extended
- **Min. Seat Height:** 23" (58.4cm)
- **Max. Seat Height:** 38" (96.5cm)
- **Overall Width:** 31.5" (80.0cm)
- **Seat Width:** 18" (45.72cm)
- **Arm Length:** 23" (58.4cm)