

Our Largest Stainless Steel Storage Cabinet



- All-Welded Construction
- Smooth Operating, Sliding Doors
- 5 Adjustable Stainless Steel Shelves
- Optional Locks for Added Security
- Adjustable Leveling Glides
- 3 Year Limited Warranty

The UMF Medical SS7816 stainless steel storage and supply cabinet is ideal for storage of valuable equipment, supplies and precision instruments. It includes five removable, adjustable stainless steel shelves, adjustable in 1/2" (1.3cm) increments, and also comes standard with glass panel doors for convenient visibility of all contents.

Optional locks are available to ensure security of valuables or restricted access items.

3 year limited warranty.

Standard Features

- Heavy duty, all-welded stainless steel construction
- Five removable, adjustable stainless steel shelves
- Full Height, stainless steel sliding doors with glass panes
- Recessed finger pulls on doors for ease of operation
- Adjustable leveling glides
- 3 year limited warranty

Specifications

- Width: 47" (119.4cm)
- Height: 84.25" (214.0cm)
- Depth: 15.875" (40.3cm)

Options

- 503 – Sliding door lock (factory installed only)
- 507 – Additional keys (x2)
- NPN – solid panel doors instead of glass panes