



Spacious Storage, Full View Doors

All-Welded Construction
Hinged, Glass Panel Doors with Lock
5 Adjustable Glass Shelves
Adjustable Leveling Glides
3 Year Limited Warranty

The UMF Medical SS7840 stainless steel storage and supply cabinet features full height stainless steel hinged doors with glass panes for convenient and secure storage as well as easy access to your precision instruments. Locking T-Handle on right side door locks both doors in the closed position. Five removable stainless steel shelves are included, adjustable in 1/2" (1.3cm) increments.

3 year limited warranty included.

Standard Features

- Heavy duty, all-welded stainless steel construction
- Five removable, adjustable stainless steel shelves
- Full Height, stainless steel sliding doors with glass panes
- Convenient T-handle door latches with lock
- Adjustable leveling glides
- 3 year limited warranty

Specifications

- Width: 35.875" (91.1cm)
- Height: 76.75" (195cm)
- Depth: 16" (40.6cm)

Options

- 507 – Additional keys (x2)
- NPN – solid panel doors instead of glass panes