

Bassinet, One Shelf with Basket and Mattress

Convenient Storage Shelf
Basket Can Be Placed Horizontally or Elevated on One End
3" (76mm) Ball Bearing Swivel Casters For Safe Transport
304 Medical Grade Stainless Steel Construction
Basket and Mattress Pad Included



Standard Features

- Heavy duty tubular stainless steel framework
- Convenient storage shelf
- 3" (7.6cm) ball bearing swivel casters
- 3 year limited warranty
- Includes basket and mattress pad
- Name card holder

Specifications

- **Height:** 35.375" (89.9cm)
- **Width:** 31.25" (79.4cm)
- **Depth:** 17" (43.2cm)

Options

- 132 - Swivel casters locking 3" (set of 4)