

WELCH ALLYN ELITE™ RETINOSCOPE



Superior Illumination in a Prestigious New Retinoscope

The Welch Allyn Elite Retinoscope's enhanced light output makes refractions easier than ever by producing one of the brightest, sharpest, and cleanest streaks available. Features such as a convenient external focusing sleeve, sealed optics, and an optional polarizing filter make the Welch Allyn Elite Retinoscope one of the most user-friendly and innovative instruments offered.

WelchAllyn®

Advancing Frontline Care™

SUPERIOR ILLUMINATION FOR UNPARALLELED PERFORMANCE

For more than 85 years, Welch Allyn diagnostic instruments have set the standards for superior performance and value. And with the introduction of the Welch Allyn Elite™ Retinoscope, we've raised the bar once again. This prestigious instrument offers a combination of outstanding features never before seen in a retinoscope. Which means you'll be able to perform examinations with greater ease and accuracy than ever before.

The improvements incorporated in the new Welch Allyn Elite Retinoscope have been achieved by making several modifications to our internal design, including aperture size changes, the use of new lenses with ant reflective coatings and internal texturing to reduce glare.

Superior Streak Quality

Enables easier and quicker refractions
– even in patients with the smallest pupils

Crossed Linear Polarizing Filter

Eliminates 99% of bothersome glare from trial lenses, helping you perform a more precise exam at the flick of a switch


External Focusing Sleeve

The continuously rotating external focusing sleeve improves comfort and manoeuvrability of the scope, no matter what size your hand is, or how you prefer to manipulate the sleeve

Dust-Free Design

Special construction keeps the optics cleaner for more effective examinations

BRIGHTNESS VS. STANDARD RETINOSCOPE



Elite Optical System

The Elite Optical System provides 2.8X the brightness and superior streak quality versus standard retinoscopes.

US Patents

5,189,449 • 5,200,772 • 5,202,710

Lamp Patent

5,359,262

Spot Retinoscopy Interchangeability

Simply insert a spot filament bulb in place of the streak filament bulb, and the Welch Allyn Elite Retinoscope is converted into a Spot Retinoscope. This allows you to perform spot retinoscopy without having to buy a second scope.

Ordering Information

18245	Welch Allyn Elite Retinoscope
71020-C	Convertible Rechargeable Power Handle
71670	3.5v Rechargeable Handle for Desk Charger
71904	Lithium Ion Power Handle with AC Charging Module - UK
71902	Lithium Ion Power Handle with AC Charging Module- Europe
08200-U	Halogen HPX™ Streak Replacement Lamp
08300-U	Halogen HPX Spot Replacement Lamp



European Customer Services Centre,
Welch Allyn Limited, Navan Business Park, Dublin Road, Navan, Co Meath, Republic of Ireland
Tel. 00-353-469067790 Fax. 00-353-469067754

Welch Allyn UK Limited, Cublington Road, Aston Abbots, Buckinghamshire, HP22 4ND
Tel. 0207-365-6780 Fax. 0207-365-9694
website : www.welchallyn.co.uk

SM2231EU RevC

WelchAllyn®

Advancing Frontline Care™