

Presenting a Convenient, Affordable Diagnostic System that Leaves Little Room for Improvement...



A Complete, Affordable Diagnostic System

The Welch Allyn 767 Integrated Diagnostic System gives you all the instruments you'll need to perform a complete diagnostic examination.



767 Wall Transformer

Handle-based rheostats and automatic on/off sensors make operation simple.



MacroView[™] Otoscope

A nearly complete view of the tympanic membrane with 30% greater magnification than a traditional otoscope.



767 Wall Aneroid Sphygmomanometer

Easy-to-read, jeweled-movement, mercury-free manometer is certified to be accurate within one percent of scale (±3 mmHg).



Pneumatic Otoscope

Rotating lens with ultra seal for pneumatic otoscopy.



PanOptic™ Ophthalmoscope

5x larger field of view of the fundus vs. standard ophthalmoscopes in an undilated eye (25° vs. 5°).

11810—6 Filters/Apertures

11820—6 Filters/Apertures including: Cobalt Blue Filter and Add-On Corneal Magnifying Lens



KleenSpec® Plus Disposable Specula Dispenser

Small/large—optimum storage and organization of specula and supplies.



3.5 V Direct Ophthalmoscope

Long-lasting, reliable performance for general and specialist examinations.

11710—5 Apertures and Red-free Filter

11720—Coaxial optical system including: 6 Apertures and 3 Filters



SureTemp® Plus Electronic Thermometer

Provides accurate oral temperature readings in about 4 to 6 seconds and rectal temperature readings in about 10 seconds.

And Plenty of Room for Everything Else.




WelchAllyn®









767 Integrated Diagnostic System Components Chart

COMPONENTS PRODUCT #	Wall Board	767 Transformer 76710	Traditional Ophthalmoscopes		PanOptic™ Ophthalmoscopes		Diagnostic MacroView™ Otoscopes	Pneumatic Otoscopes	767 Aneroid	KleenSpec® Dispenser (small)	KleenSpec® Dispenser (large)	SureTemp® Plus Thermometer
			11710	11720	11810	11820	23810	20200	7670-01	52401-F	52101-F	01690-300
76791-MX	•	•	•				•		•	•		
76791-1MPX	•				•		•		•	•		
76791-2MPX	•					•			•	•	-14 [7]	W •
76792-M	•	•	•				•		•	•	4-7-1	1011
76792-1MP		•							•	•	1,1413	
76792-2MP	•	•				•	•		•	•		
76793-MX	•	•					•		•	•		•
76794-M	•	•					•		•	•		
76797-X								•	•		•	•
76797-2PX	•	• -				•		•	•		•	•
76798	•	a • E						•	•		•	
76798-2P	•							•	•		•	

767 Integrated Diagnostic System for Spot Vital Signs® or Spot Vital Signs® LXi

Wall Board allows easy mounting of Spot or Spot LXi devices to Integrated System. Order Spot or Spot LXi separately.

COMPONENTS	Wall Board	767 Transformer	Traditional Ophthalmoscopes		PanOptic™ Ophthalmoscopes		Diagnostic MacroView™ Otoscopes	Pneumatic Otoscopes	767 Aneroid	KleenSpec® Dispenser (small)	KleenSpec® Dispenser (large)	SureTemp® Plus Thermometer
PRODUCT #		76710	11710	11720	11810	11820	23810	20200	7670-01	52401	52101	01690-300
76794-MSPOT	•	•					•			•		
76794-MLXI	•	•		•			•			•		
76794-2MPSPOT	•	•					•			•		
76794-2MPLXI	•	•				•	•			•		