



RIGID ENDOSCOPY

WelchAllyn[®]

Advancing Frontline Care[™]

Quality and Value in

Welch Allyn offers a wide range of rigid endoscopy products for anoscopy and sigmoidoscopy, complemented by versatile fibre optic illumination systems for optimum visualisation.

The sigmoidoscope range includes both reusable and disposable instruments and a variety of sizes to accommodate differing requirements. Welch Allyn fibre optic anoscopes represent the finest in design and craftsmanship for optimal performance.



Welch Allyn Disposable Sigmoidoscopes

One of our most popular designs, the Welch Allyn Disposable Sigmoidoscope may be used in conjunction with either a 6 V illumination system or traditional Welch Allyn 3.5 V power handles.

- > Highly visible graduations for greater accuracy and control
- > Glare-free frosted coating ensures maximum light output at distal tip
- > Water-activated lubricating gel to facilitate easy insertion
- > Individually wrapped for safety and convenience
- > Fibre-optic illumination system provides a cool 360 degree ring of light
- > Hinged window with neoprene seal to avoid fogging during insufflation

53130	KleenSpec® Disposable Sigmoidoscope with Obturator (case of 100), (add 'L' suffix for lubricated version)
73500	3.5 V Halogen Handle Adapter
71020-A	3.5 V Rechargeable Handle
73210	Light Handle with 0.9 m Cord
36019	Fibre Optic Light Head
73322	Transformer with 1.5 m Cord (EUR)

Insufflation Systems

An insufflation system can be used with Welch Allyn sigmoidoscopes to ease insertion of the instrument. Hydrophobic filters are available to prevent contamination of the insufflation bulb. They are attached using a Luer lock fitting.



30200	Insufflation Bulb, complete
30215	Luer Lock for use with hydrophobic filters (pack of 10)
31210	Hydrophobic filters (pack of 10)

Rigid Endoscopy

Standard and Paediatric/Stricture Sigmoidoscopes

All Welch Allyn reusable sigmoidoscopes feature:

- > Autoclavable stainless steel for convenience and durability
- > Hinged window with a neoprene seal to avoid fogging during insufflation
- > Graduated tube with removable tip for fibre cleaning (metal only)



Standard Sigmoidoscopes

32810	19 mm x 15 cm Standard Sigmoidoscope, Speculum only with Obturator
32820	19 mm x 25 cm Standard Sigmoidoscope, Speculum only with Obturator
33830	19 mm x 30 cm Standard Sigmoidoscope, Speculum only with Obturator
32830	19 mm x 35 cm Standard Sigmoidoscope, Speculum only with Obturator
33220	23 mm x 25 cm Standard Sigmoidoscope, Speculum only with Obturator



Fibre Optic Sigmoidoscopes

32010	11 mm x 10 cm Fibre Optic/Stricture Speculum only with Obturator
32020	11 mm x 25 cm Fibre Optic/Stricture Speculum only with Obturator
32410	15 mm x 15 cm Fibre Optic/Paediatric Speculum only with Obturator
32420	15 mm x 25 cm Paediatric/Stricture Speculum only with Obturator

Welch Allyn Anoscopes

A complete line of high-performance and reusable anoscopes for safe and efficient examinations.



Disposable Anoscope

53110 KleenSpec® Disposable Anoscope with Obturator (box of 25)

Disposable Anoscope Accessories

73210	Light Handle with 0.9 m Cord
73222	Light Handle (73210) plus Transformer (73322)
36019	Fibre Optic Light Head
73322	Transformer with 1.5 m Cord (EUR)



Rotating Anoscope

Anoscope Only 68 mm Speculum Length

38900 Rotating Speculum Anoscope with Obturator, 20 mm inside diameter, 22 mm outside diameter




Fibre Optic Anoscope

Speculum Only with Obturator

37019	19 mm x 7 cm Fibre Optic Anoscope, Speculum Only with Obturator
37023	23 mm x 7 cm Fibre Optic Anoscope, Speculum Only with Obturator
37027	27 mm x 7 cm Fibre Optic Anoscope, Speculum Only with Obturator

Welch Allyn reusable anoscopes are also available in Long-Speculum, Operating and Standard configurations. Please contact us for further information on sizes and product codes.

Flexible Illumination Systems for Optimal Visualisation



Sigmoidoscope/Anoscope Lighting Systems

Powerful illumination systems.

- > Halogen light provides consistent, long-lasting illumination (73500, 73550 only)
- > Fibre optics provide a cool distal ring of light with no reflections or obstructions (36019, 38700 only)
- > Sigmoidoscopes/anoscopes can be used with a choice of illumination systems—a 6 V system with optional on/off switch and rheostat control, or with Welch Allyn's 3.5 V/2.8 V or 2.5 V hand-held power sources (with adaptors 73500 and 73550).
- > 73500 adapter fits most Welch Allyn 3.5 V power handles
- > 73550 adapter fits most Welch Allyn 2.5 V power handles (except PocketScope™)

Sigmoidoscope/Anoscope Lighting Systems and Accessories

36019	Fibre Optic Light Head for Disposable Sigmoidoscopes and Anoscopes	73550	2.5 V Halogen Handle Adapter
38300	2.5 V Vacuum Light Carrier	30200	Insufflation Bulb, Complete
38700	Fibre Optic Light Carrier for use with Standard Anoscopes and Light Adapter	30130	30 cm Suction Tube (Standard)
73210	Light Handle with 0.9 m cord	30140	40 cm Suction Tube
73308	Optional On/Off Switch	34000	Rectal Probe
73322	Transformer with 1.5 m cord (EUR)	34100	Rectal Hook
73400	Optional Rheostat Control	00200-U	2.5 V Vacuum Lamp (small) (for 38300)
73500	3.5 V Halogen Handle Adapter	03100-U	3.5 V Halogen Lamp (for 73500)
		03400-U	2.5 V Halogen Lamp (for 73550)
		07800-U	6 V Halogen Lamp (for 73210)

For more information, contact us at +353 46 906 7790 today!

Welch Allyn European Customer Service, Navan Business Park,
Dublin Road, Navan, County Meath, Ireland
Tel: +353 46 906 7790 Fax: +353 46 906 7754
www.welchallyn.com

Welch Allyn PTY. LTD, S. Africa
394 Surrey Ave, Ferndale, Randburg 2194 South Africa
P O Box 809, Fontainebleau, South Africa 2032
Tel: 0 11 777 7555 Fax: 0 11 777 7556
www.welchallyn.com

WelchAllyn®

Advancing Frontline Care™